
SINCE 1966
MADE IN ITALY

NebulAIR⁺

by **FLAEM**[®]

GB INSTRUCTIONS FOR USE MANUAL

NO BRUKSANVISNING

Assembly diagram - Monteringskjema

NebulAIR⁺

Mod. NEBULAIR

AEROSOL THERAPY APPARATUS

We are pleased you have purchased our product and we thank you for your trust in us. We aim at fully satisfying our customers by offering them state-of-the-art products for the treatment of respiratory diseases.

Read these instructions carefully and retain them for future reference. Only use the accessory as described in this manual. This is a medical device to nebulise and administer medication prescribed or recommended by your doctor.

Please note that the full range of

Flaem products is visible on the website www.flaem.it

STANDARD EQUIPMENT INCLUDES:

Ⓐ Aerosol apparatus (main unit)

- A1 - On/Off switch
- A2 - Air intake
- A3 - Air filter
- A4 - Nebuliser holder
(internal/external)
- A5 - Handle for transport
- A6 - Cable storage compartment
- A7 - Power cable

Ⓒ Accessories

- C1 - RF7 Dual Speed Plus Nebuliser
 - C1.1 - Lower part
 - C1.2 - Nozzle
 - C1.3 - Upper part
 - C1.4 - Speed selector with valve
- C2 - Mouthpiece with valve
 - C2.1 - Exhalation valve
- C3 - Non-invasive nosepiece
- C4 - Paediatric SoftTouch mask
- C5 - Adult SoftTouch mask
- C6 - Manual nebulisation control
- C7 - Filter replacement kit

Ⓑ Connection tube (main unit / nebuliser)

⚠ IMPORTANT WARNINGS

- Before using the product for the first time, and periodically during its lifetime, check the integrity of the device structure and of the power cable to make sure there is no damage. In the event of damage, do not plug in the cable and immediately take the product to an authorised FLAEM service centre or to your trusted dealer.
- The expected medical life of the accessories is minimum 1 year when used 4 times a day.
- Children and people who are not self-sufficient must use the device under the close supervision of an adult who has read this manual.
- Some parts of the device are small enough to be swallowed by children; therefore, keep the device out of the reach of children.
- Do not use the supplied tubing and cables for anything other than their intended use. These parts could cause a strangling hazard: pay close attention to children and persons with particular difficulties as they are often unable to accurately evaluate danger.
- The apparatus is unsuitable for use in presence of flammable anaesthetic mixture with air, oxygen or nitrous oxide.
- Always keep the power supply cable away from hot surfaces.
- Keep the power cable away from animals (for example, rodents) which could damage the insulation.
- Do not handle the device with wet hands. Do not use the device in damp environments (for example, while taking a bath or shower). Do not immerse the device in water; in the event of immersion immediately disconnect the plug. Do not remove or touch the

immersed device; unplug the power cable first. Immediately bring the device to an authorised FLAEM service centre or to your trusted dealer.

- Use the device only in dust-free conditions, otherwise treatment could be compromised.
- The device casing is not waterproof. Do not wash the device under running water or by immersion and keep it safe from being sprayed by water or other liquids.
- Do not expose the device to particularly extreme temperatures.
- Do not place the device near sources of heat, in direct sunlight or in excessively hot rooms.
- Do not obstruct or put objects into the filter or its related housing in the device.
- Never obstruct the air vents located on both sides of the device.
- Always use it on a rigid surface that is clear of obstacles.
- Make sure there is no material obstructing the air vents before each use.
- Do not put any objects in the air vents.
- Repairs must be carried out only by FLAEM authorised personnel. Unauthorised repairs void the warranty and may be hazardous for the user.
- The average expected duration for the compressor series are: F400: 400 hours, F700: 700 hours, F1000: 1000 hours, F2000: 2000 hours.
- WARNING: Do not modify this device without authorisation from the manufacturer.
- The Manufacturer, the Vendor and the Importer shall be held responsible for safety, reliability and performance only if: a) the device is used in compliance with the instructions for use b) the wiring where the device is being used is in compliance with safety regulations and current laws.
- Interactions: the materials used in contact with medication have been tested with a vast range of medications. However, in view of the variety and continuous evolution of pharmaceuticals, interactions cannot be ruled out. We recommend using the medication as soon as possible once it has been opened and preventing prolonged exposure in the nebuliser cup.
- The manufacturer must be contacted to report problems and/or unexpected events related to operation.

OPERATING INSTRUCTIONS

Before each utilisation, wash your hands accurately and clean the device as described in the paragraph "CLEANING, SANITISATION, DISINFECTION AND STERILISATION".

The nebuliser cup and accessories are for personal use only to avoid potential risks of transmitting infection.

During application it is advisable to suitably protect yourself from any dripping.

1. Plug the power cable (A7) into mains socket that is compatible with the device voltage. The position of the socket must be such that the device can be easily unplugged from the mains network.
2. Insert the nozzle (C1.2) in the upper part (C1.3) pressing as shown by the 2 arrows in the "Connection diagram" in point C1. Insert the Speed selector with valve (C1.4) in the upper part (C1.3) as shown in the "Connection diagram" in point C1. Pour the medication prescribed by the doctor into the lower part (C1.1). Close the nebuliser by turning the upper part (C1.3) clockwise.
3. Connect accessories as indicated in the "Connection diagram" on the cover.
4. Sit comfortably holding the nebuliser in your hand, place the mouthpiece onto your mouth or alternatively use the nose piece or mask. Should you opt for the mask acces-

sory, place it on your face as shown in the picture (with or without using the elasticated strap).

5. Start the device by means of switch (A1) and breathe deeply in and out. After inhaling, we recommend holding your breath for an instant to allow the inhaled drops of aerosol to deposit. Then exhale slowly.
6. Upon completing application, switch off the device and unplug it.


ATTENTION: If after the therapy session an evident deposit of moisture forms within the pipe (B), detach the pipe from the nebuliser and dry it with the very ventilation from the compressor; this operation prevents possible blooms of mould inside the pipe.

USE METHODS OF THE RF7 DUAL SPEED PLUS NEBULISER WITH SPEED SELECTION AND VALVULAR SYSTEM


It is professional, quick, recommended to administer all types of medications, including the costlier ones, even in patients with chronic diseases. Thanks to the geometry of the RF7 Dual Speed Plus nebuliser cup internal ducts, we have obtained a recommended, active particle size to treat even the lower respiratory tract.

For faster inhalation therapy, position the valve speed selector (C1.4) pressing on the word MAX with your finger.

For more effective inhalation therapy, position the valve speed selector (C1.4) pressing on the opposite side of the word Max with your finger. In this case, you have optimal assumption of the medication, minimising leaks into the surrounding environment, thanks to the valvular system with which the nebuliser cup, mouthpiece and mask are equipped.


SOFTTOUCH MASKS

The **SoftTouch** masks have an external edge made in **soft, biocompatible material** that guarantees optimal adherence to the face and are also equipped with the **innovative Dispersion Limiting device**. These distinctive features allow increased sedimentation of the medication within the patient's body, also **limiting its dispersion**.


During inhalation, the tab that works as a **Dispersion Limiting device** bends toward the inside of the mask.


During exhalation, the tab that works as a **Dispersion Limiting device** bends toward the outside of the mask.

USE OF NEBULIZATION MANUAL CONTROL

To achieve continuous nebulisation action you should not use the manual nebulisation control (C6), especially in the case of children or persons with reduced physical, sensory, or mental capabilities. The manual nebulisation control is useful for limiting dispersion of the medication in the surrounding environment.

To start nebulizing close with a finger the hole of the nebulizer manual control (C6) and breathe in gently; we recommend to hold your breath for a moment so that the inhaled aerosol droplets can be deposited,

meanwhile, to disable nebulizing, remove your finger from the hole of the nebulizer manual control to avoid the waste of drug, optimizing its acquisition. Then exhale slowly.

CLEANING, SANITISATION, DISINFECTION, STERILISATION

Switch off the device before any cleaning procedure and unplug the power cable from the socket.

DEVICE AND TUBING EXTERIOR

Use only a damp cloth with antibacterial soap (non-abrasive and with no solvents of any sort).

ACCESSORIES

Open the nebuliser by turning the upper part (C1.3) anticlockwise, remove the nozzle (C1.2) from the upper part (C1.3) by pressing as shown by the 2 arrows in the "Connection diagram" in point C1.

Then proceed according to the following instructions.

CLEANING AT HOME - SANITISATION AND DISINFECTION

SANITISATION

Before and after each use, sanitise the nebuliser cup and the accessories, choosing one of the methods described below.

method A: Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent.

method B: Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** in the dishwasher with a hot cycle.

method C: Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** by immersing them in a solution of 50% water and 50% white vinegar, then rinse thoroughly under potable hot water (approximately 40°C).

If you want to also perform the cleaning for DISINFECTION, jump to the DISINFECTION paragraph.

After having sanitised the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, a hair dryer).

DISINFECTION

After sanitising the nebuliser cup and the accessories, disinfect them choosing one of the methods described below.

method A: Accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** can be disinfected.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

Implementation:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.

- Remove the disinfected parts and rinse abundantly with warm potable water.
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

method B: Sanitise the accessories **C1.1-C1.2-C1.3-C1.4-C2-C3** by boiling them in water for 20 minutes; use demineralised or distilled water to prevent calcium deposits.

method C: Sanitise the accessories **C1.1-C1.2-C1.3-C1.4-C2-C3** with a hot steam steriliser for baby-bottle (not the microwave type). Perform the process faithfully following the instructions of the steriliser. To ensure that the disinfection is effective, choose a steriliser with an operating cycle of at least 6 minutes.

After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, from a hair dryer).

CLEANING IN A CLINICAL OR HOSPITAL SETTING - DISINFECTION AND STERILISATION

Before disinfection or sterilisation, sanitise the nebuliser cup and the accessories, choosing one of the methods described below.

method A: Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** under potable hot water (approximately 40°C) with a gentle, non abrasive dish detergent.

method B: Sanitise accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** in the dishwasher with a hot cycle.

DISINFECTION

Accessories **C1.1-C1.2-C1.3-C1.4-C2-C3-C4-C5-C6** can be disinfected.

The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

Implementation:

- Fill a container big enough to hold all of the parts to disinfect with a solution of potable water and disinfectant, according to the proportions indicated on the packaging of the disinfectant.
- Completely immerse each part in the solution, taking care to avoid the formation of air bubbles on the parts. Leave the parts immersed for the amount of time indicated on the packaging of the disinfectant associated with the concentration chosen for the solution.
- Remove the disinfected parts and rinse abundantly with warm potable water.
- Dispose of the solution following the instructions provided by the disinfectant manufacturer.

If you want to also perform the STERILISATION, jump to the STERILISATION paragraph.

After having disinfected the accessories, shake them vigorously and place them on a paper towel. Alternatively, dry them with a jet of hot air (for example, from a hair dryer).

STERILISATION

Accessories **C1 (nebulizer assembled) - C2-C3-C4-C5** can be sterilised.

Device: Fractionated vacuum overpressure steam steriliser in accordance with EN 13060.

Implementation: Wrap every single part to be treated with a sterile barrier system or packaging in accordance with Norm EN 11607. Place the packed components in the steam steriliser. Run the sterilisation cycle according to the operating instructions of the device by selecting a temperature of 134°C and a time of 10 minutes first.


Storage: Store the sterilised parts as per the instructions for use of either the sterile barrier system or packaging.

The sterilisation procedure is validated in its conformity to ISO 17665-1.

AIR FILTERING

The device is quipped with an extraction filter (A3) that must be replaced when it is dirty or changes colour. Do not wash or reuse the same filter. Regular filter replacement is necessary to help ensure proper compressor performance. The filter must be checked regularly.

Contact your retailer or authorised service centre for replacement filters.

To replace the filter:

Insert a slotted screwdriver between the edge of the filter and the body. Lift the filter and extract it by rotating and pulling up. The filter was designed to always remain secured in its seat. **Use only with original Flaem accessories.**

SYMBOLS

Type Approved

www.tuv.com
ID 2000000000

TUV approval ref. EN 60601-1 3rd Ed. - EN 60601-1-11

EC Marking medical ref. Dir 93/42
EEC and subsequent updates

Class II device

Attention: check the instructions
for use

Operating switch off

Operating switch on

In compliance with: European
Standard EN 10993-1 "Biological
Evaluation of medical devices"
and European Directive 93/42/EEC
"Medical Devices."
Phthalate-free. In compliance with:
Reg. (EC) no. 1907/2006

Device serial number

Manufacturer

Type BF applied part

Alternating current

Risk: electrocution.
Consequence: Death.
Using the device while taking a bath
or shower is prohibited

Degree of protection of the casing:
IP21.
(Protected against solid foreign objects larger
than 12mm. Protected against access with a fin-
ger; Protected against vertically dripping water).

SPARE PARTS

Description

- RF7 Dual Speed Plus Nebuliser cup set, comprising:
 - nebuliser cup, mouthpiece with valve and nosepiece
- 2 m tubing connection
- Masks set, comprising:
 - adult mask, paediatric mask and elastic band
- Air filter replacement kit (2pcs)

TECHNICAL SPECIFICATIONS

Mod. NEBULAIR

Power source:	230V ~ 50Hz 210VA
Max pressure:	3.5 ± 0.5 bar
Compressor air delivery:	approx 14 l/min
Noise (at 1 m):	approx 55 dB (A)
Operation:	Continuous
Operating Conditions:	
Temperature:	min 10°C; max 40°C
Air humidity:	min 10%; max 95%
Atmospheric pressure:	min 69KPa; max 106KPa
Storage conditions:	
Temperature:	min -25°C; max 70°C
Air humidity:	min 10%; max 95%
Atmospheric pressure:	min 69KPa; max 106KPa

Safety approvals:

Dimensions:	20 (W) x 30 (D) x 10 (H) cm
Weight:	2.400 Kg

APPLIED PARTS

Type BF applied parts are: patient accessories (C2, C3, C4, C5)

RF7 Dual Speed Plus Nebuliser

Minimum medication capacity:	2 ml
Maximum medication capacity:	8 ml
Operating pressure (with neb.):	approx. 1.30 bar

	C1.4 speed selector	
	in Max pos	in Min pos with valve
(1) Delivery:	approx. 0.65 ml/min	approx. 0.29 ml/min
(2) MMAD:	2.38 µm	2.21 µm
(2) Breathable fraction < 5 µm (FPF):	81.5%	84.2%

(1) data recorded according to internal Flaem procedure I29-P07.5.

(2) In vitro characterization certified by TÜV Rheinland LGA Products GmbH - Germany in compliance with European Standard for Aerosol therapy apparatuses EN 13544-1. Further details are available on request.

DEVICE DISPOSAL

In compliance with the Directive 2012/19/EC, the symbol printed on the device shows that the device to be disposed of is considered waste and must therefore be an item of "differentiated collection". Consequently, the user must take it (or have it taken) to the designated collection sites provided by the local authorities, or turn it in to the dealer when purchasing an equivalent new device. Differentiated waste collection and the subsequent treatment, recycling and disposal procedures promote the production of devices made with recycled materials and limit the negative effects on the environment and on health caused by potential improper waste management. Unauthorised disposal of the product by the user shall result in the application of administrative penalties, as per article 50 and the subsequent revisions of Italian Legislative Decree no. 22/1997.

ELECTROMAGNETIC COMPATIBILITY

This device was designed to satisfy the currently required requisites for electromagnetic compatibility (EN 60601-1-2:2007). Electro-medical devices require particular care during installation and use relative to EMC requirements. Users are therefore requested to install and/or use these devices following the manufacturer's specifications. There is a risk of potential electromagnetic interference with other devices, in particular with other analysis and treatment devices. RF mobile or portable radio and telecommunications devices (mobile telephones or wireless connections) can interfere with the functioning of electro-medical devices. For further information visit our website www.flaemnuova.it. Flaem reserves the right to make technical and functional modifications to the product with no prior warning.

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© 2017 FLAEM NUOVA®
FLAEM NUOVA S.p.A.